

Store PDF files to the PACS




PacsFiler

PacsFiler™ provides any Windows user the option to file important PDF documents in the PACS. By using PacsFiler these documents can be linked to existing patient information and archived in the safety of the PACS

Easy to use

PacsFiler™ is easy to use and does not require any special training or DICOM knowledge. The workflow is simple. Just launch the application and in three steps, your document is permanently linked to the patient study, archived in your PACS.



Step 1. Select a Patient Study

Selection of the desired patient and study starts the filing session. Each PDF document will be saved in the PACS as a single image.

New patients can also be created at this time.

Step 2. Select the document to be filed

Use the import button to select the document which should be filed for the selected patient.

Step 3. Store

Just click store and the document will be filed in the PACS.

That is all !

Features

- MWL, Patient Root and Study root queries¹
- New patients can be added
- View loaded PDF in default PDF viewer
- Available for Windows XP Professional, Vista and Windows 7²

1) The DICOM conformance statement for PacsFiler is available and provides essential information for DICOM specialists and others with an interest in the technical details of this product.

2) Windows 2000 support is available on request.

We highly recommend PacsPrinter (another Medima product) for printing/filing documents to the PACS that are not available in the PDF format

©2010 Medima

Interested? Go to www.medima.nl and download a free demo version!