

The document scanning solution for your PACS




PacsScanner™ scans a paper document and stores the pages as a series of images in your PACS. These documents can be anything: insurance cards, lab reports, handwritten notes, etc. In this way all patient information is conveniently available at every PACS workstation for instant retrieval.

Easy to use

PacsScanner™ is easy to use and does not require any DICOM knowledge. The workflow is simple.

Select a patient study

The patient list is populated with studies from the PACS or Worklist. Searching is made easy by selection of studies through common demographic filters.




Name	Patient ID	BirthDate	Sex
김희중	2008-3	18000101	M
Buc^Jérôme	SCSFREN	19010101	M
ヤマダ太郎=山田太郎=やまだたろう	H32EXAMPLE	19730621	F
Äneas^Rüdiger	SCSGERM	20010227	M
Hong^Gildong=洪^吉洞=홍^길동	I2EXAMPLE	19511025	M
Wang^XiaoDong=王^小东	X2EXAMPLE	19170303	M
Wang^XiaoDong=王^小東	X1EXAMPLE	20010606	F
Διονυσιος	SCSGREEK	19531025	M
Локсембур	SCSRUSS	19661115	M
שרון דבורה	SCSHBRW	19581011	M
قباڻي^المنار	SCSARAB	19430505	M
ヤマダ太郎=山田太郎=やまだたろう	H32EXAMPLE	19470430	F

Study ID	Description	Date	Accession nu...
SCSGREEK	CHEST	20081213	4981038

Interested? Go to www.medima.nl and download a free demo version!

- Scan the document and finally press the Store button

All the most common scanning options are present; images can be flipped with full support for black and white and color images. More pages can be included in the same image series. Pictures and/or movies can be also be loaded from the hard disk.



- Scan paper documents** Any paper document can be imported from the scanner
- Acquire images from web-cameras** Pictures from any connected web-cam can be acquired and archived
- Import images and movie** Picture and movies can be loaded from the hard disk and archived
- No extra workstation required** PacsScanner™ works on standard Windows PCs networked to the PACS system
- Simple user interface** Anyone can scan documents without any DICOM or PACS knowledge
- Document manipulation** PacsScanner™ allows flipping of documents and images and selection of resolution before transmission
- Support for international charactersets** PacsScanner™ recognizes and displays all the fonts: arabic, chinese, cyrillic, etc.

©2010 Medima

- Interested? Go to www.medima.nl and download a free demo version!