

Print anything to the PACS



PacsPrinter


PacsPrinter™ provides any Windows application with PACS connectivity. Just by pressing the *Print* button on any Windows application your document will be sent to the PACS.

Easy to use

PacsPrinter™ is easy to use and does not require any DICOM knowledge. The workflow is simple.

Select the correct printer

After selecting the *Print* option of the application you are working with, you can choose PacsPrinter as destination.



Select Patient and Study Demo Mode

Options Help

Search Add Patient

Name: Patientid: Birthdate: Sex:

Accession Nr.: Studydate:

Name	Patient ID	BirthDate	Sex
Buc^Jérôme	SCSFREN	19000101	M
Hong^Gildong=洪^吉洞=홍^길동	I2EXAMPLE	19531025	M
Wang^XiaoDong=王^小东	X2EXAMPLE	20010606	F
Wang^XiaoDong=王^小東	X1EXAMPLE	20010606	F
Yamada^Tarou=山田^太郎=やまだ^たろう	H31EXAMPLE	19740430	M
Διοσσοδς	SCSGREEK	19030203	M
Люксембур	SCSRUSS	19661115	M
قباني^اليزار	SCSARAB	19430505	M
IMAGPATIENT1209^1209	000-00-1209	19171231	M
IMAGPATIENT567^567	000-00-0567	19430000	M
שרון^דבורה	SCSHBRW	19340221	M
ヤマダ^太郎=山田^太郎=やまだ^たろう	H32EXAMPLE	19740430	M
조희준	2008-3	18000101	M

IMAGPATIENT1209^1209
 Patientid : 000-00-1209 Birthdate : 12-31-1917 Sex : M

Study ID	Description	Date	Accession nu...
test	test	20091202	GMRC-123

Study ID: Study Description: Ref.Physician:

Series number: Series Description: Operator name:

medima logo and buttons: Search, Clear, Don't print, Continue

Select a patient study

Selection of the desired patient and study completes the printing session, documents will be saved into the PACS as a series of images. New patients can be created. No additional software has to be used: you can access PacsPrinter™ from *any* Windows application.

PacsPrinter 2.0 Feature Matrix

Feature	PacsPrinter Deluxe	PacsPrinter Premium	PacsPrinter Pro
Archive from any Windows program	•	•	•
Allow new patients to be added	•	•	•
Query Option			
Modality Worklist Query	•	•	•
Patientroot Query	•	•	•
Studyroot Query	•	•	•
Dicom output			
Secondary Capture	•		•
Encapsulated PDF		•	•
Supported Operating Systems			
XP	•	•	•
Vista	•	•	•
Windows 7	•	•	•